

ALWAYS USE MULTIPLE REAGENTS TO AVOID FALSE POSITIVES - THERE ARE MORE SUBSTANCES THAN COLORS.

QUICKLY SEARCH RESULTS USING PROTESTKIT.EU/RESULTS WEB APP.

Results will vary depending on sample adulteration, quantity, lighting, temperature, exposure to air, storage, as well as reagent quality and degradation. Positive or negative result does not indicate if a substance is safe. No substance is 100% safe.



- Canabidiol (CBD) - main test for THC / CBD, helps rule out synthetic cannabinoids
Erlisch - main test for all opiates, e.g. lysergamides and tryptamines but also 5-HTP and more, secondary test for cocaine
Hofmann - main test for lysergamides and tryptamines, more precise but slower than Erlisch
Froehde - main test for most powders, secondary test for ketamine and amphetamines, long shelf life
Liebermann - main test for most powders, secondary test for cocaine (detects levamisole)
Marquis - main test for most powders, secondary test for ketamine, methadone and psychedelics
Mecke - main test for most powders, secondary test for cocaine, ketamine and methadone
Morris - main test for cocaine and ketamine, can tell apart some novel dissociatives
Rohdoffe - main test for primary amines, e.g. MDMA, amphetamine or PMMA
Simons - main test for secondary amines, e.g. MDMA, methy, methadone, methylone or PMMA
Zimmermann / Beam - main test for cathinones and benzodiazepines in general, detects also CBD
Fentanyl Test Strip to detect fentanyl, even if it is mixed in another substance

Table comparing reagents with columns for various substances and rows for different tests, indicating results with symbols like +, -, or +/-. Includes a legend for symbols.

REAGENT CALIBRATION

- ERLISCH: 1. Light, 2. Red
HOFMANN: 1. Blue, 2. Red
FROEHDE: 1. Purple, 2. Red
LIEBERMANN: 1. Purple, 2. Red
MARQUIS: 1. Purple, 2. Red
MECKE: 1. Purple, 2. Red
ROHDOFFE: 1. Purple, 2. Red

CANNABIS AB REAGENT RESULTS

- THC > CBD: Your sample is THC dominant
CBD > THC: Your sample is CBD dominant, if it has strong effects it is adulterated!
Different result?: Risk of false negative result by using more samples...

Table showing calibration results for ERLISCH, HOFMANN, FROEHDE, LIEBERMANN, MARQUIS, MECKE, and ROHDOFFE reagents.

ERLISCH REAGENT RESULTS

Table of ERLISCH reagent results for various substances including opiates, amphetamines, ketamine, and others.

MORRIS REAGENT RESULTS

Table of MORRIS reagent results for cocaine, ketamine, and other substances.

FENTANYL TEST STRIP RESULTS



Testing method #1: If testing a pill, crush one in a plastic baggie and dump out the powder...
Testing method #2: Dissolve 100 mg of a sample in 200 ml of water...

FROEHDE REAGENT RESULTS

Table of FROEHDE reagent results for a wide range of substances.

HOFMANN REAGENT RESULTS

Table of HOFMANN reagent results for various substances.

MORRIS REAGENT RESULTS

Table of MORRIS reagent results for cocaine and ketamine.

SIMONS REAGENT RESULTS

Table of SIMONS reagent results for secondary amines.

ZIMMERMANN REAGENT RESULTS

Table of ZIMMERMANN reagent results for cathinones and benzodiazepines.

LIEBERMANN REAGENT RESULTS

Table of LIEBERMANN reagent results for cocaine and other substances.

MARQUIS REAGENT RESULTS

Table of MARQUIS reagent results for ketamine, methadone, and psychedelics.

MORRIS REAGENT RESULTS

Table of MORRIS reagent results for cocaine and ketamine.

SIMONS REAGENT RESULTS

Table of SIMONS reagent results for secondary amines.

ZIMMERMANN REAGENT RESULTS

Table of ZIMMERMANN reagent results for cathinones and benzodiazepines.

MARQUIS REAGENT RESULTS

Table of MARQUIS reagent results for ketamine, methadone, and psychedelics.

MORRIS REAGENT RESULTS

Table of MORRIS reagent results for cocaine and ketamine.

SIMONS REAGENT RESULTS

Table of SIMONS reagent results for secondary amines.

ZIMMERMANN REAGENT RESULTS

Table of ZIMMERMANN reagent results for cathinones and benzodiazepines.

MECKE REAGENT RESULTS

Table of MECKE reagent results for cocaine, ketamine, and methadone.

MORRIS REAGENT RESULTS

Table of MORRIS reagent results for cocaine and ketamine.

SIMONS REAGENT RESULTS

Table of SIMONS reagent results for secondary amines.

ZIMMERMANN REAGENT RESULTS

Table of ZIMMERMANN reagent results for cathinones and benzodiazepines.